

FDA Syringe Recall History

- FDA recalls for dispensable products:
 - Cardinal Health Monoject syringes – [link 1](#) , [link 2](#)
 - Medline syringes – [link](#)
 - Jiangsu Shenli Medical Production Co. Ltd – [link](#)
 - Sol-Millennium Medical, Inc – [link](#)
 - Fresenius Medical Care – [link](#)
 - B. Braun Medical, Inc - [link](#)

BD Alaris System v12.3 Syringe Pump Compatibility

System v9.33 Syringe Compatibility

- BD Plastipak:
 - 1, 3, 5, 10, 20, 30, and 50/60 mL
- Covidien Monoject:
 - 3, 6, 12, 20, 35, and 60 mL
- Terumo:
 - 3, 5, 10, 20, 30, and 50/60 mL
- Astra Zeneca:
 - 50 mL
- IVAC
 - 50 mL

System 12.3 Syringe Compatibility

- BD:
 - 1, 3, 5, 10, 20, 30, and 50 mL
- Covidien Monoject:
 - 12, 35, and 60 mL
- Terumo:
 - 12 mL

BD Alaris System v12.3 PCA Pump Compatibility

System v9.33 PCA Syringe Compatibility

- BD Plastipak:
 - 20, 30, and 50/60 mL
- Covidien Monoject:
 - 20, 35, and 60 mL
- Terumo:
 - 20, 30, and 50/60 mL
- IMS:
 - 30 mL

System 12.3 PCA Syringe Compatibility


- BD:
 - 20, 30, and 50 mL
- Covidien Monoject:
 - 35 and 60 mL

503b PCA Syringe Information

- Currently available information (as of 4/5/2024):
 - QuVa:
 - Per representative, utilize Covidien Monoject PCA syringes
 - Leither:
 - Per representative, utilize Covidien Monoject PCA syringes
 - SCA:
 - Pending response from representatives

Alternative Syringe Types

- Alternative syringe types:
 - [Aerogen](#) continuous nebulization syringe
 - Utilized by respiratory therapy for continuous nebulization administration
 - Per Aerogen, this is a BD syringe rebranded as Aerogen and replaced blue plunger
 - Compatible with BD Alaris
 - [EquaShield](#) CTSD syringes
 - Utilized for hazardous and/or chemotherapy administration
 - Per manufacturer, the syringe is fully manufactured in house
 - Not validated nor recognized by BD for syringe pump use
 - [Medline](#) syringes
 - Generic luer lock syringe
 - Per Medline, the syringe is a Covidien Monoject syringe rebranded as generic Medline
 - Not validated nor recognized by BD for syringe pump use



ICU-Medical Medfusion 4000 Syringe Pump Syringe Compatibility

- BD Plastipak:
 - 1, 3, 5, 10, 20, 30, 50, and 60 mL luer-lock
 - 1 mL tuberculin
- Covidien Monoject:
 - 3, 6, 12, 20, 35, and 60 mL luer-lock
 - 1 mL tuberculin
- Terumo:
 - 3, 5, 10, 20, 30, and 50/60 mL luer-lock
 - 1 mL tuberculin
- B. Braun:
 - 5, 10, 20, and 50 mL luer-lock